



SKR054-ET77A

EMERGENCY TROLLEY

Register code: SKM-A



SAIKANG Medical

Healthier
your life

SKR054-ET77A

FEATURES AND HIGHLIGHTS

Technical parameters

- External size (LxWxH) 750×475×950mm

Technical configuration

- | | |
|------------------------|------|
| ▪ Defibrillator shelf | 1pc |
| ▪ Sliding Side Shelf | 1pc |
| ▪ Sharp Container | 1pc |
| ▪ Dust Baskets | 2pc |
| ▪ IV pole | 1pc |
| ▪ CPR Board | 1pc |
| ▪ Oxygen Bottle Holder | 1pc |
| ▪ Power Outlet | 1pc |
| ▪ Central Lock | 1set |
| ▪ File Box | 1pc |



Material

Steel and aluminum structure inside, ABS and antibacterial materials outside, easy to clean.



Extendable Worktop

Laterally extendable worktop




Table-Board

With handles on both sides, surface with transparent soft glass




Drawer

Inner with partitions can be well organized



IV Pole

Movable IV pole can exchange position with defibrillator shelf

**Defibrillator Shelf**

The plate surface is made of virgin ABS material. The position of the IV pole and the defibrillator shelf can be exchanged for easy use.

**Power Outlet**

The power socket is movable and detachable, Easy to use, Located above the CPR board to convenient the use of oxygen cylinders

**Oxygen Bottle Holde**

Inserted oxygen bottle holder, belts fixed on trolley main body

**Trash Can**

Two ABS trash cans with lid

**Debris Basket**

With two debris baskets on both sides. The big one is for two sharp containers and the small one is for sanitizer

**Wheels**

TPR tire no worn out after running 30KM, save anti-winding hard shell, United forming without bolts
Pass Dynamic test: Bearing 120kg run 30KM and pass obstacles 500times





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service

WORLDWIDE BRANCHES

South America Branch

Saikang Medical Technology, Inc
W3400 Corrientes, Argentina
TEL: +55 (11) 2758-0438
Rua Universal, Hollywood,
São Bernardo do Campo,
Brazil
TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd
Mezzanine Floor, Taj Tower,
upper hill, nairobi, Kenya
TEL: +254714084806

South East Asia Branch

Saikang Medical Technology Ltd
Sala Thammasop Subdistrict
Thawi
Watthana District
Bangkok Thailand
TEL: +66 642741046

UK Branch

Uk Healthward International Limited
DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +31 (0) 2021111106

Russia Branch

LLC 'SAIKANG'
454004, Russian Federation,
Chelyabinsk,
Universitetskaya
Naberezhnaya, Russia
TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies (Mid East) LLC
Musaffah
P.O Box 53200, Abu Dhabi,
UAE
Tel: +971 585710652

SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
 +86 18021231901 · 24/7
 export@saikangmedical.com
 www.saikangmedical.com



Company video



General catalogue