

**mimsal**

# LUXIFLEX LED II

**Examination** and **advanced exploration** luminaire.

## Version

Luxiflex LED II  
**15.000 lux at 50 cm**



## High efficiency

The powerful LED light source offers exceptional light quality and a lifespan of more than 60,000 hours.

## Flexibility

Flexible, height-adjustable and highly precise arm that allows the luminaire to be adjusted to any position.

## Design & handling

Designed with clean and soft shapes which facilitate hygiene and cleanliness in medical environments.

Its compact head and the switch located at the base of the lamp make this model an ideal luminaire for the most demanding examinations procedures.

## Mounting options

The standard version includes AH table clamp.

As for accessories, we offer a wide selection of supports such as wall supports, wall extension arms, trolley stands and rail brackets.



# LUXIFLEX LED II

| Technical Details          | (REF. LX2370) |
|----------------------------|---------------|
| Light source               | LED 8W        |
| Number of LEDs             | 3             |
| Useful life                | 60.000 hours  |
| Illumination at 50 cm      | 15.000 lux    |
| Color temperature          | 5.500 °K      |
| Color rendering            | > 97%         |
| R <sub>9</sub> Color index | > 90%         |
| Total length               | 935 mm        |
| Color                      | White         |

## ACCESSORIES



Wall B Support  
REF. B0101800



Ah Table Clamp  
REF. B010190



Rail Plus Bracket  
REF. 82026



Extension Arm  
REF. B0101S0



Trolley stand  
4 kg REF. 09655  
8,8 kg REF. 09651