




Veterinary scale

6858

- ▶ Quick and easy operation via four function keys
- ▶ 4 adjustable feet for levelling
- ▶ Units in kg, lb, oz, lb.oz
- ▶ Energy-saving auto-off function



Features

- ▶ Quick and easy operation via four function keys
- ▶ Easy-to-read LC display, digit height 25 mm
- ▶ Tare function: documents placed on top can be tared
- ▶ Hold function for better readability
- ▶ Aluminium display unit
- ▶ Large stainless steel weighing platform with separate grey rubber mat
- ▶ 4 adjustable feet for levelling
- ▶ Easy to transport thanks to handle and two castors
- ▶ Units in kg, lb, oz, lb.oz
- ▶ Energy-saving auto-off function
- ▶ Mains operation DC 5V mains adapter
- ▶ Optional battery operation with four 1.5V AA batteries (batteries not included in the scope of delivery)

Products

Order number Variation	Approvable	Division g	Maximum load kg	Platform size mm
6858.01.002 Veterinary scale	No	50	200	900 x 600