SPHYGMOMANOMETERS AND STETHOSCOPES STETHOSCOPES LATEX-FREE



HEINE GAMMA® STETHOSCOPES

- : Solid chest piece and ideal weight avoid stray noise from manipulation and ensure excellent resonance. Hygienic matt chrome finish.
- : Chest piece with non-chill rims for the comfort of the patients.
- : Ergonomic, chrome-plated binaurals. Individually adjustable.
- : High-quality Y-piece tubing for perfect acoustics and sound reproduction. Length 56 cm.
- : Latex-free materials.

HEINE GAMMA® 3.1 Pulse Stethoscope

For fine acoustic performance.

- : Stethoscope with solid flat chest piece. The chest piece can be easily inserted under the sphygmomanometer cuff.
- : High-quality membrane, 44 mm dia.
- : Single-tube construction with excellent sound transmission.
- : Includes ear olives (large/soft) for sound insulation.
- : Optional available: Ear olives, (small/hard) (M-000.09.946).

GAMMA 3.1

HEINE GAMMA® 3.2 Acoustic Stethoscope

For excellent acoustic performance.

Technical data as GAMMA 3.1, but with a double chest piece.

- **: Combined double chest piece** with two different frequency ranges to choose from: A high-quality membrane on one side and a conical acoustic head on the other side.
- : High-quality membrane, 44 mm diameter for excellent acoustic performance and sound transmission.
- : Cone diameter 30mm.

GAMMA 3.2

M-000.09.942

M-000.09.941

HEINE GAMMA® 3.3 Acoustic Stethoscope

For the examination of children.

Technical data as GAMMA 3.2.

: Reduced contact area of 31 mm makes GAMMA 3.3 ideal for examining children.

GAMMA 3.3

M-000.09.943



HEINE GAMMA® C3 Cardio Stethoscope

Specially-developed for cardiology.

- : Exceptional acoustic performance, particularly in the cardiac frequency range.
- : **Combined double chest-piece** with two different frequency ranges. Instant selection of a double membrane on one side or a conical acoustic head on the other.
- : High-performance double membrane, 45 mm diameter gives ideal sound transmission and perfect skin contact for exceptional sound reproduction.
- : **Two-in-one tube system** for enhanced performance. Very supple, yet kink-resistant with exceptional acoustic reproduction.
- : Includes ear olives (large/soft) and one pair of replacements (small/hard) for sound insulation.

GAMMAC3

