

SCRUB-UP UNITS



Made entirely of AISI 304 stainless steel, 18/10, and Scotch Brite brushed finish. Welding works are polished, smoothed and refined with the original used stainless steel sheet.

Anti-splash -shaped basin with anti-splash guard on 3 sides. It is sound -dampened at the bottom. Stainless steel drain siphon directly welded to the basin not to have holes, and prepared to be connected to the drain trap.

Rear box with flowing top to avoid dust deposits, that contains complete plumbing and electric system to deliver water to the faucet.

Scrub-up units are available as follows:

- ▶ one station – two stations – three stations (length : 600 – 1200 – 1800 mm);
- ▶ elbow, electric photocell, battery photocell, knee – actuated faucets and mixers, complete with thermostatic pre - mixer for water temperature adjustment;
- ▶ available accessories: detergent and disinfectant dispensers, brush dispensers and hand hot air dryers.

Complete supply of installation accessories and brackets to fix the scrub-up unit to the wall in addition to all electric and hydraulic connections.

