


## Blood-pressure cuffs and accessories


### Latex-free blood-pressure cuffs

Cuffs for GAMMA G7/G5/GP/XXL LF sphygmomanometers.



- ⌘ Latex-free.
- ⌘ Fit the most widely-used sphygmomanometers.
- ⌘ High-quality velcro cuffs made of wear-resistant material. Physiologically tested.
- ⌘ Double-sewn.
- ⌘ Washable.
- ⌘ Index mark for exact positioning.
- ⌘ Clear indication of size and control mark for arm circumference.
- ⌘ Extra-thick, highly-flexible tube. For a long life.
- ⌘ Proven reliability. Absolutely airtight system.

	GAMMA G7/G5/GP			Single-tube cuffs (cpl.)		Single-tube bladders		Outer cuffs	
	Size	Circumference/cm		Size	Size	Size	Size		
Small child	6x28	8–13	<b>M-000.09.700</b>	5x15	<b>M-000.09.661</b>	6x28	<b>M-000.09.059</b>		
Child	10x34	13–20	<b>M-000.09.702</b>	8x20	<b>M-000.09.663</b>	10x34	<b>M-000.09.060</b>		
Small adult	14x53	20–29	<b>M-000.09.704</b>	11x23	<b>M-000.09.664</b>	14x53	<b>M-000.09.061</b>		
Adult	14x58	29–41	<b>M-000.09.619</b>	12x28	<b>M-000.09.665</b>	14x58	<b>M-000.09.055</b>		
Large adult	20x86	>41	<b>M-000.09.615</b>	18x52	<b>M-000.09.667</b>	20x86	<b>M-000.09.057</b>		

	GAMMA XXL LF			Twin-tube cuffs (cpl.)		Twin-tube bladders		Outer cuffs	
	Size	Circumference/cm		Size	Size	Size	Size		
Small child	6x28	8–13	<b>M-000.09.701</b>	5x15	<b>M-000.09.671</b>	6x28	<b>M-000.09.059</b>		
Child	10x34	13–20	<b>M-000.09.703</b>	8x20	<b>M-000.09.673</b>	10x34	<b>M-000.09.060</b>		
small adult	14x53	20–29	<b>M-000.09.705</b>	11x23	<b>M-000.09.674</b>	14x53	<b>M-000.09.061</b>		
Adult	14x58	29–41	<b>M-000.09.629</b>	12x28	<b>M-000.09.675</b>	14x58	<b>M-000.09.055</b>		
Large adult	20x86	>41	<b>M-000.09.625</b>	18x52	<b>M-000.09.677</b>	20x86	<b>M-000.09.057</b>		