

HbA1c results in minutes

Informed conversations at the point of care



A1CNow+®

The A1CNow+® system from Chek Diagnostics provides healthcare professionals with a fast and easy way of getting HbA1c results—allowing informed conversations about how patients are managing their diabetes in minutes, not days.

For more information, visit www.a1cnow.com.

Key Features

- National Glycohemoglobin Standardization Program certified
- CLIA waived
- Reimbursable

To place your order, call 1-877-870-5610.

 **CHEK** Diagnostics
Innovating Point-of-Care

HbA1c results in minutes

Informed conversations at the point of care



A1CNow+®

More efficient than the lab

- Results in 5 minutes
- Small (5µL) blood sample
- Portable – use in multiple exam rooms
- Easy to use – minimal training required
- Reimbursable



A1CNow+® system
(10-count test kit)
3024



A1CNow+® system
(20-count test kit)
3021

To place your order, call 1-877-870-5610.

a1cnow.com / chekdiagnostics.com

© 2014 PTS, Inc. A1CNow is a trademark of Polymer Technology Systems, Inc. LIT-001641 Rev. 0 03/14

 **CHEK** Diagnostics

Innovating Point-of-Care