


Lighted Articulating Ear Curette™



The **Lighted Articulating Ear Curette** combines the two greatest innovations in ear curettage, visualization with light and magnification, and articulation for the ability to work behind an occlusion in the ear canal. The easy to maneuver tip articulates up to 90° to assist in the removal of cerumen from the ear canal and foreign bodies from the ear and nose.

► PROBLEM:

Cerumen plugs and foreign bodies in the ear are often difficult to remove. Practitioners lack the ability to visualize or work from behind the occlusion to remove it.

► SOLUTION:

The Lighted Articulating Ear Curette brilliantly illuminates the ear canal and provides an articulating tip that allows the practitioner to get behind the occlusion.

► RESULTS:

Direct visualization and an articulating ear curette eliminates the blind, top down sweeping motion that can abrade the ear canal improving safety and efficiency.

PRODUCT OFFERING

#2511 – Lighted Articulating VersaScoop® (3mm tip)

Includes 25 Articulating Ear Curette Tips,
1 Light Source and 1 Magnification Lens.

#2200 – Extra Light Source

#2265 – Pack of 5 Magnification Lenses

POTENTIAL REIMBURSEMENT (U.S. ONLY)

CPT Code / Description

69210 - For potential reimbursement of impacted cerumen.

69200 - For potential reimbursement of foreign body removal from the ear.

30300 - For potential reimbursement of foreign body removal from the nose.



©2014 Bionix Medical Technologies. Patent Pending.



Obelis s.a • Boulevard Général Wahis 53 • 1030 Brussels, BELGIUM
Tel: +(32) 2.732.59.54 • Fax: +(32) 2.732.60.03 • E-Mail: mail@obelis.net



Bionix Development Corp. • 5154 Enterprise Blvd. Toledo, OH, U.S.A.
Phone: 800-551-7096 • Fax: 800-455-5678 • www.Bionix.com

GOING BEYOND ►

Phone: 800.551.6810
www.BionixMed.com

