

HEINE® mPack mini

Compact and lightweight power source for LoupeLight 2 / MicroLight 2



- EXCLUSIVELY AT HEINE**
- Operating time typ. 9 hours
 - Stepless light intensity control
 - Compact and lightweight

The mPack mini is only compatible with LoupeLight 2 and MicroLight 2.

- **Freedom from cables** that tie the user to desk, wall or static power supplies.
- **Operating time typ. 9 hours.** In continuous operation at full power typ. 4 hours.
- **Charging time typ. 4 hours.**
- **Charge status indicator.** Display of remaining battery power with three LEDs. The flashing of the last LED indicates that immediate charging is required. During charging the LEDs light up one after the other.
- **Stepless light intensity control.** Optimum brightness setting between 100% and 3% reduces reflexes.
- **USB cord** with E4-USB Medical approved plug-in power supply.
- **The integrated galvanic separation module** protects the battery and the user against the consequences of defective power supplies, thus preventing the possibility of the full mains voltage flowing directly to the housing. We recommend the use of a medical approved plug-in power supply.
- **Simple cleaning.** Closed housing, no seams at the sides, few edges.
- **Compact and lightweight.** With a weight of just 81g and dimensions of 44 x 23 x 101mm the power source can be easily stowed.
- **mPack mini Belt Clip.** Pouch to easily secure the power source to belt or clothing.
- **Li-ion technology.** No „memory effect“, without worrying about charge level or battery condition.



[01]



[02]



[03]

mPack mini	6V
mPack mini with Li-ion rechargeable battery and E4-USB plug-in transformer	X-007.99.650
mPack mini with Li-ion rechargeable battery, without E4-USB plug-in transformer	X-007.99.649
mPack mini Belt Clip [01]	X-000.99.007
USB cord [02] with E4-USB plug-in transformer [03]	X-000.99.303
E4-USB Medical approved plug-in power supply for USB cord [03]	X-000.99.305