

Pulse Oximeter

PO1000

Tabletop Pulse Oximeter

- High resolution 4.3 inch color LCD display
- Compact and portable for uninterrupted monitoring
- AC power supply with built-in lithium battery
- Adjustable audible and visual alarms
- Storage and review up to 36 hours
- USB data upload
- Optional for printing



PO2000

Tabletop Pulse Oximeter

- Color LCD display and screen rotation
- Automatic power-off function for power saving
- Spo2 and Pulse rate measurement
- Three selectable working modes
- Suitable for adult, pediatric and neonatal patients
- Data can be transferred to PC for storage review and printing
- Long battery up to 20h working on 4 X AA alkaline or Ni-MH rechargeable batteries
- Powerful data storage capacity
- Audible and Visible alarm capability



PO80

Fingertip Pulse Oximeter

Features :

- Dual Color OLED display with four directions
- Compact and power-saving
- Battery capability indication
- Spo2, PR, Spo2 bar graph and pulse waveform display
- Power supply: two standard AAA batteries
- Automatic switch-off function for power-saving
- Dimension: 58mm*36mm*32mm

Standard Configuration: SPO2, PR
Optional Configuration: Menu function

Color:

Blue, Grey, Purple, Green, Pink



PO81

Fingertip Pulse Oximeter

Features :

- Four directions and Six display modes
- Spo2, PR, Spo2 bar graph and pulse waveform display
- Power supply: two 1.5V AAA batteries
- Operating menu for function setting
- Power consumption: below 30mAh
- Dimension: 58mm*35mm*30mm

Standard Configuration: SPO2, PR
Optional Configuration: USB, Menu function

Color:

Black, White, Golden



PO82

Fingertip Pulse Oximeter

Features :

- Dual color OLED display with various display modes
- Low power consumption and indication
- Automatic switch-off no later than 8 seconds after no signal
- Optional function: Automatic power-on after 5 seconds when the product detect the signal
- Optional function: Gravity sensor function. Four different display modes will change along with the movement of hands
- Optional function: PI value which is used to reflect the situation of pulsation flow.

Color:

Blue, Green, White, Orange

