

## Wall- / Table-Transformer HEINE® EN 50

### Modular Power Supply System/charger



⚡ Modular Power Supply System

⚡ Soft Start Technology

⚡ Switches off automatically



The HEINE Modular Power Supply System offers unsurpassed flexibility and innovative new features. The EN50 powers OMEGA500 and ML4 LED HeadLight. Furthermore, all other HEINE 6V products are compatible. The EN50 Charger, the mPack and HC50 L Headband Mounted Rheostat can be combined in every conceivable configuration.

- ⚡ **Modular, compact design** functions as the base for Tabletop or Wall Mounted use and as charging base for the mPack mobile charging system.
- ⚡ **Upgradeable.** mPack Charging base by simple removal of the control module.
- ⚡ **Soft Start Technology.** A gradual introduction of current reduces bulb “shock” and extends bulb life.
- ⚡ **100 – 240V Compatible Transformer** incl. 4 plugs for world wide compatibility.
- ⚡ **Registered Design.**

EN50 Charger	6V
EN50 Transformer table/wall unit with Rheostat Control Module (incl. 90° Angled Adaptor, Extension Cord Cinch-Cinch 2m, Instrument support) [01]	X-095.17.301
Rheostat Control Module for EN50 [02]	X-095.17.305
EN50 Base unit (without Rheostat Control Module) [03]	X-095.17.300
Plug-in 90° Angled Adaptor (for EN50/mPack without instrument support) [04]	C-000.33.509
Instrument support for Headband	C-000.33.511