

Smokerlyzer[®] breath carbon monoxide (CO) monitors



CE
0086



Helping people to stop smoking, one breath at a time.

www.bedfont.com



NEW

iCO™ Smokerlyzer®

Personal use

The world's first Smokerlyzer® breath CO monitor that connects to your smartphone



Features and benefits

- Your personal CO monitor
- The only Smokerlyzer® you can fit in your pocket
- Connects to the Smokerlyzer® app on your phone or tablet
- Send your results to your doctor, smoking cessation advisor, family member or friend
- Monitors your progress, encouraging you to quit smoking.

INSTANTLY see your personal CO results

Technical specification

Concentration range:	0-100ppm
Detection principle:	Electrochemical Sensor
Accuracy (repeatability of reading):	≤± 5 %
Hydrogen cross-sensitivity:	≤6%
Operating life:	Approx. 200 tests/ 3 years (whichever comes first)
Sensor sensitivity:	1ppm
MOQ:	Please contact Bedfont®
Dimensions:	Approx. 28mm diameter x 115mm long
Weight:	Approx. 35g
System requirement:	iPad, iPhone and iPod require iOS 6 and above. Android phones and tablets require Android OS Version 4.0
Warranty:	6 months

- Ideal for:**
- Personal use¹¹
 - Health professionals
 - Pharmaceutical companies
 - App developers
 - Insurance companies



Contact Bedfont® or one of our worldwide **Smokerlyzer®** distributors for a free demonstration.

www.bedfont.com
Tel: +44 (0)1622 851122
Email: ask@bedfont.com

A full list of our worldwide distributors can be found at
www.bedfont.com/distributors/

Our family, innovating health, for yours.

Bedfont® Scientific Ltd
Station Road, Harrietsham, Maidstone,
Kent, ME17 1JA England
Tel: +44 (0)1622 851122 Fax: +44 (0)1622 854860
Email: ask@bedfont.com
www.bedfont.com



© Bedfont® Scientific Limited 2017

Issue 17- March 2017, Part No: MKT261
Bedfont® Scientific Limited reserves the right to change or update this literature without prior notice.
Registered in: England and Wales, Registered No: 1289798