



# Pulse Oximeter

For Home Use

Cherish life,keep healthy

# JPD-500B

## Pulse Oximeter



- ❖ Integrated with SPO<sub>2</sub> probe and processing display module.
- ❖ Small in volume, light in weight and convenient in carrying.
- ❖ Operation of the product is simple, low power consumption.
- ❖ SPO<sub>2</sub> and PR display,
- ❖ Pulse rate value display, bar graph display.
- ❖ Low-voltage indication: low-voltage indicator appears before working abnormally which is due to low-voltage.
- ❖ Automatically power off function: when the device is under the state of measuring interface.
- ❖ Various color of cover can be selected.
- ❖ Auto shut off :16s
- ❖ LED display
- ❖ SPO<sub>2</sub> : 35~100%
- ❖ Pulse Rate : 25~250BPM
- ❖ Resolution: 1% for SPO<sub>2</sub>, 1BPM for pulse rate
- ❖ Accuracy: ±1% (70%~100%),  
unspecified (<70%) for SPO<sub>2</sub>,  
±1BPM or ±1% (select larger) for pulse rate
- ❖ Power: 1.5V (AAA size) alkaline batteries × 2
- ❖ Size: 58.0×31.6×32.9mm



**Shenzhen Jumper Medical Equipment Co.,Ltd.**

---

Address: Jumper Building, Fuyong First Industrial Zone, Xintian Road  
Fuyong Street, Bao'an District, Shenzhen, 518103 P.R.China

Tel.: +86-755-26696279

Fax: +86-755-26852025

Website: [www.jumper-medical.com](http://www.jumper-medical.com)

E-mail: [info@jumper-medical.com](mailto:info@jumper-medical.com)