

MIMLED

Compact LED Surgery lamp.

Available in two versions: **MIMLED 1000** and **MIMLED 600**.

HIGH EFFICIENCY

MIMLED 1000 100.000 lux @ 1 meter

MIMLED 600 60.000 lux @ 1 meter

HIGH COLOUR RENDERING

The high colour rendering $Ra > 94\%$, describes the exceptional quality of white light, allows tissues to be distinguished with utmost precision.

LIGHT INTENSITY DIMMING

Adjustable light intensity integrated in the handgrip and on the head which can be regulated using the panel control.

ENERGY SAVING LIGHT

Low power consumption.

33W MIMLED 1000

25W MIMLED 600

AVAILABLE IN 3 OPTIONS

- Trolley Stand.
- Wall version.
- Ceiling version (up to 3.10 meter height).



MIMLED 1000



MIMLED 600



MIMLED 1000

Technical Details	
Lux at 1 meter	100.000 lux
Lamp lifetime	50.000 hours
Colour temperature	4.500 °K
Colour rendering	≥ 94%
Light field diameter	195 mm
Power supply	100/230 V - 50/60 Hz
Consumption	33 W

MIMLED 600


Technical Details	
Lux at 1 meter	60.000 lux
Lamp lifetime	50.000 hours
Colour temperature	4.500 °K
Colour rendering	≥ 94%
Light field diameter	195 mm
Power supply	100/230 V - 50/60 Hz
Consumption	25 W



*Adjustable light intensity
integrated in the grip*



Ceiling version
Cod. ML1000C
ML600C



Wall version
Cod. ML1000W
ML600W



Trolley Stand version
Cod. ML1000FL
ML600FL

