

## EU TYPE EXAMINATION CERTIFICATE

**Certificate No: 2163-PPE-1902**

Respiratory protective devices, filtering half masks to protect against particles manufactured by

**Jiangmen Yanyang Trading Co., Ltd.**

No.1, 4th Floor, Building 2, No.18 Xinyi Road, Jianghai District, Jiangmen City, Guangdong  
Province, CHINA

are tested and evaluated according to

**EN 149:2001 + A1:2009 Respiratory Protective Devices -  
Filtering Half Masks to Protect Against Particles -  
Requirements, Testing, Marking**

Based on the type examination conducted with the evaluation of test reports, technical file  
according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 5, it is approved  
that the product meets the requirements of the regulation.

**Product Definition**

Single shift use particle filtering half mask for protection against solid and liquid aerosols, is a  
folding type, 5 layers, without valve, nose clip, head straps and internal sponge.

**Model: YY1029**

**Classification: FFP3 NR**

For more details, refer technical evaluation report provided to the manufacturer, dated  
05.01.2021 and number 2163-KKD-1902

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as  
shown below, on the Category III product models given above, with;

- Issuing an appropriate EU Declaration of Conformity according to **Personal Protective Equipment Regulation (EU) 2016/425 Annex 9.**
- Ongoing successful performance in fulfilment of the requirements set out in **Personal Protective Equipment Regulation (EU) 2016/425** and harmonised standards, ensured by assessments based on **Annex 7 (Module C2) or Annex 8 (Module D)** of the regulation.

This certificate is initially issued on **09/01/2021** and will be valid for 5 years, if there is no change  
in the relevant harmonised standard affecting the essential health and safety requirements.



  
**Suat KAÇMAZ**  
UNIVERSAL CERTIFICATION  
Director

