



Fitness scale

7850

- Measurement of muscle, fat and water components, BMI and caloric needs per day
- Printout of measurement values
- Simple self-explanatory menu guide on the glass surface
- Clothing weight value seasonally adjustable



SOEHNLE PROFESSIONAL	
24-03-2010 09:50	
Height:	160cm
Age:	26
Sex:	female
Weight with clothing:	55.4kg
Weight without clothing:	52.9kg
BMI: Gewicht ok	20.7

Features

The fitness scale is ideal for use in fitness studios, pharmacies, nutrition counselling and wellness hotels. With its body analysis function, in addition to weight and body mass index, the muscle, fat and water content can also be measured precisely and the current energy metabolism can be determined. The measurement is taken by hand, so that the user can stand on the scales comfortably and save time while wearing shoes. Compared to currently known analysis scales, the scale is particularly easy to use and does not require complicated and detailed instructions. With the integrated printer, the weight and analysis values including individual evaluation can be conveniently printed out, documented and easily tracked over a longer period of time.

- ▶ Measurement of muscle, fat and water components in the body
- ▶ Impedance measurement via the hands, no need to take off the shoes
- ▶ Large LC display
- ▶ Indication of time in standby mode
- ▶ Simple self-explanatory menu guide on the glass surface
- ▶ No limit on measurement for over- and underweight people
- ▶ Clothing weight value seasonally adjustable
- ▶ Automatic data printout of the analysis measurements
- ▶ Optional: Company information on printout individually adjustable prior to delivery (see accessories)
- ▶ Extra large platform: 320 x 400 mm
- ▶ Switchable to kg / st / lb units
- ▶ Mains operation 100 - 240 V, 50-60 Hz

Available languages:

German, English, French, Spanish, Italian, Greek, Turkish, Hungarian, Dutch, Czech, Danish, Finnish, Swedish, Norwegian, Slovak, Polish, Russian

Products

Order number Variation	Approvable	Division g	Maximum load kg	Platform size mm
7850.30.001 Fitness scale	No	100	180	365 x 490

Accessories

Order number	Description
Variation	
2500.01.008	Paper for thermal printer (10 rolls, 58 mm)
5751.01.007	Customization of the printout
2557.25.001	UK adapter