

➤ **HEINE mini 3000®**

Compact pocket instruments – Maximum quality with minimum dimensions.

NEW! HEINE miniNT® charger set



➤ **HEINE mini 3000®**

The complete range of high-quality diagnostic instruments in a modern and ergonomic design.



MINI TO THE MAX. The handy “mini” series is the name of the compact range of HEINE pocket instruments with the high-quality features of HEINE’s “full-size” instruments. Characteristics include an innovative and compact design, coupled with high-quality, full-range features for the maximum in diagnostic reliability. The series also features automatic power shut-off when returned to the pocket, together with superb stability and manufacturing.

CONVENIENTLY-SIZED
AND COMPACT.

Bigger doesn’t always mean better. In the case of the mini3000 series, the unique ergonomics and the high-quality features represent the crucial characteristics of superlative quality. The result is exceptional handling in all investigation situations.

COMPLETE AND PRECISE.

Our most important goal when developing the mini3000 series was optimised technical sophistication and enhanced performance. We’ve delivered this through innovative features and improved XHL Xenon Halogen Lens bulb technology for better diagnostic quality and reliability.

MODERN,
ERGONOMIC DESIGN.

Quality meets design – and with the mini3000 series, the phrase couldn’t be more true. Alongside a modern design, special plastics ensure supreme sturdiness yet low weight. The “mini” series also satisfies every demand in terms of functionality and performance.



HEINE mini 3000®
BATTERY HANDLE

HEINE mini 3000®
OPHTHALMOSCOPE

HEINE mini 3000®
FOCALUX

HEINE mini 3000®
TONGUE-BLADE HOLDER

Compact battery handle
suitable for all mini 3000
instruments.



Innovative, compact battery
handle in a high-quality,
modern design.

- ✦ **High-quality handle:**
Chrome-finish upper section
/refined plastic. Shockproof,
sturdy, non-slip.
- ✦ **Oval, compact, stylish shape.**
Lies ergonomically in the hand.
- ✦ **Attachment clip with**
integrated on/off switch.
Secure. Switches off
automatically when returned
to the pocket.
- ✦ **Screw-on end cap.** For easy
battery change.
- ✦ **Threaded connector made**
of metal. Prevents wear.
- ✦ **Complete range of mini 3000**
instruments. Versatile range.
- ✦ **Replaceable batteries.**
Size AA.

Ophthalmoscope in a unique,
compact design with
outstanding optics for even
the smallest dimensions.
Perfect when on the move
and doing home visits.

- ✦ **Enhanced XHL Xenon**
Halogen Technology.
Concentrated white light for
perfect illumination and a
brilliant image.
- ✦ **5 different apertures.**
A full-function instrument incl.
fixation star.
- ✦ **Range of lenses:**
-20 to +20 D.

Focusing lamp for the
examination of the anterior
chamber in focal light and for
indirect ophthalmoscopy.

- ✦ **Small and light.** Particularly
ideal for home visits.
- ✦ **Light spot can be focused**
on the patient's pupil.
Perfect fundus illumination
even with non-dilated pupil.
- ✦ **Enhanced XHL Xenon**
Halogen Technology for
more light. Very bright,
concentrated white light with
high light intensity for a bright
fundus image.
- ✦ **+3 D lens attachment.**
Aids accommodation and
magnifies image.

Tongue-blade holder with blade
ejector for HEINE disposable
tongue blades.

- ✦ **Enhanced XHL Xenon Halogen**
Technology. Concentrated white
light. Concentrated illumination of
oral cavity and pharynx.
- ✦ **Finger-tip pressure releases the**
used blade. No need to touch
the blade after use.
- ✦ **Solid metal construction.**
Durable and tough.
- ✦ **Only for HEINE Disposable**
Tongue-Blades.

HEINE mini 3000®
COMBI LAMP



HEINE mini 3000®
LARYNGEAL MIRROR



HEINE mini 3000®
**FIBER OPTIC (F.O.)
OTOSCOPE**



HEINE mini 3000®
OTOSCOPE



HEINE mini 3000®
DERMATOSCOPE



Combination pocket lamp and tongue blade. Durable, multi-purpose instrument.

- ✦ **Many applications.** Can be used as a general illuminator and as an illuminated tongue-blade holder for wooden blades and HEINE disposable blades.
- ✦ **Enhanced XHL Xenon Halogen Lens-lamp with 100 % more light.** Bright, focused light.

Sturdy, all-purpose laryngeal mirror for optimum illumination of oral cavity and pharynx.

- ✦ **Enhanced XHL Xenon Halogen Technology.** Concentrated illumination of oral cavity and pharynx.
- ✦ **Solid metal construction.** Durable and tough.

Otoscope with Fiber optic (F.O.) illumination.

- ✦ **Distal Fiber optic (F.O.) illumination with enhanced XHL Xenon Halogen Technology producing 100 % more light.** Homogenous, reflex-free illumination of ear canal and tympanum.
- ✦ **Viewing window with 3x magnification** and optimised casing surface for razor-sharp images and minimal reflection. Swivelling viewing window built into instrument.
- ✦ **Swivelling viewing window built into instrument.** Useful for instrumentation, cannot be mislaid.
- ✦ **Insufflation port.** For pneumatic testing of tympanic mobility. Please order connector separately.

Otoscope with direct illumination.

- ✦ **Enhanced XHL Xenon Halogen Technology with 100 % more light.** Very bright, concentrated white light for perfect illumination.
- ✦ **Viewing window with 3x magnification** and optimised casing surface for razor-sharp images and minimal reflection.
- ✦ **Swivelling viewing window built into instrument.** Useful for instrumentation, cannot be mislaid.

Dermatoscope with top-quality optics. 10x magnification with focusing optics.

- ✦ **10x magnification.** Large, high-resolution image.
- ✦ **Enhanced XHL Xenon Halogen Technology for more light.** Very bright, concentrated white light for perfect illumination.
- ✦ **Detailed compendium with photographic atlas included.** Simplifies the differential diagnosis of pigmented lesions.

HEINE mini3000®
CLIP LAMP



HEINE mini-c®
CLIP LAMP



HEINE mini-c®
CLIP LAMP WITH
EAR LIGHT



HEINE miniNT®
CHARGER- AND ADAPTOR SET



Sturdy, long-life pocket diagnostic lights in a compact, modern design with chrome head. Bright, concentrated light thanks to novel XHL Xenon Halogen lens bulb.

- ✦ **Compact, oval, attractive light.** Lies ergonomically in the hand.
- ✦ **Enhanced XHL Xenon Halogen lens bulb with 100 % more light.** Bright, focused light.

Sturdy, long-life pocket diagnostic light in a compact, modern design with chrome head. mini-c clip lamp only fits mini-c handles.

- ✦ **Compact, oval, attractive light.** Lies ergonomically in the hand.
- ✦ **Enhanced XHL Xenon Halogen lens bulb with 100 % more light.** Bright, focused light.
- ✦ **High-quality handle:** Chrome-finish upper section /refined plastic. Shockproof, sturdy, non-slip.
- ✦ **Switches off automatically when replaced in the pocket.**
- ✦ **Replaceable batteries.** Size AAA.

Clip Lamp with **Ear Light for otoplasty.** The Ear Light is used to ensure that hearing aids fit correctly.

miniNT charger set 2,5V

Consists of the miniNT charger, 2 miniNT rechargeable batteries and 2 bottom inserts for the mini3000 rechargeable handles. With this set, existing mini3000 handles can be updated to the 2.5V rechargeable system. Optional adaptor set for mini2000 available.

miniNT charger:
Width 112 mm, Height 50 mm, Depth 52 mm, Weight 270 g.

- ✦ **Innovative, compact design. Practical and spacesaving.**
- ✦ **Two independant charging compartments.** Complete functionality for one or two handles at a time. Two handles can be charged at a time.
- ✦ **Fast charge feature.** 4 hours to recharge a mini NiMH 2.5V rechargeable battery.
- ✦ **Automatic charging.** Independent of the state of charge of the rechargeable battery. Guarantees extended battery life.
- ✦ **Can not be overcharged.** Charger switches off automatically when the battery is fully-charged.
- ✦ **Charge status display.** Light pulses from the charging compartment confirm the charging mode.

Please note: Suitable only for mini 2.5V NiMH rechargeable battery.

➤ **HEINE mini 3000®**

NEW: mini NT CHARGER SET FOR mini 3000.



[HEINE mini 3000®]

2.5 V




mini3000 F.O. Otoscope with mini3000 battery handle and 1 set (4 pcs.) reusable tips and 5 of both 2.5 and 4 mm Ø disposable tips	D-001.70.120
mini3000 F.O. Otoscope with mini3000 battery handle and 5 of both 2.5 and 4 mm Ø disposable tips	D-001.70.110
mini3000 F.O. Otoscope with bulb, without handle and without accessories	D-001.70.106
Insufflation port connector	D-000.80.101
XHL Xenon Halogen spare bulb	X-001.88.105



mini3000 Otoscope with mini3000 battery handle, 1 set (4 pcs.) reusable tips and 5 of both 2.5 and 4 mm Ø disposable tips	D-001.70.220
mini3000 Otoscope with mini3000 battery handle and 5 of both 2.5 and 4 mm Ø disposable tips	D-001.70.210
mini3000 Otoscope with bulb, without handle and without accessories	D-001.70.206
XHL Xenon Halogen spare bulb	X-001.88.110



mini3000 Laryngeal Mirror with plane mirror 20 mm Ø, without handle	D-001.77.101
XHL Xenon Halogen spare bulb	X-001.88.057



mini3000 Tongue-Blade Holder complete with mini3000 battery handle and 5 disposable blades	D-001.74.118
Tongue-Blade Holder head with bulb, without handle	D-001.74.100
XHL Xenon Halogen spare bulb	X-001.88.037



mini-c Clip Lamp in blister pack with 2 AAA batteries and mini-c handle	D-001.73.109
mini-c Clip Lamp with Ear Light in blister pack with 2 AAA batteries	D-001.73.103
EarLight for mini-c Clip Lamp 6 pieces in blister pack	D-000.73.105
XHL Xenon Halogen spare bulb	X-001.88.108



mini3000 Clip Lamp complete with mini3000 battery handle	D-001.73.131
mini3000 clip lamp head for connection to mini3000 battery handle	D-001.73.130
XHL Xenon Halogen spare bulb	X-001.88.107



mini3000 Combi Lamp on mini3000 battery handle with 5 disposable blades	D-001.76.120
mini3000 combi lamp with bulb, without handle	D-001.76.101
XHL Xenon Halogen spare bulb	X-001.88.107



mini3000 Ophthalmoscope with mini3000 battery handle	D-001.71.120
mini3000 Ophthalmoscope without handle	D-001.71.105
XHL Xenon Halogen spare bulb	X-001.88.106



mini3000 Focalux with +3 D lens with mini3000 battery handle	D-001.72.141
mini3000 Focalux without lens attachment with mini3000 battery handle	D-001.72.139
mini3000 Focalux without lens, without mini3000 battery handle	D-001.72.100
Swivel lens attachment +3 D only	C-000.14.109
XHL Xenon Halogen spare bulb	X-001.88.042



mini3000 Dermatoscope with contact plate, mini3000 battery handle and 10 ml dermatoscopy oil	D-001.78.107
mini3000 Dermatoscope with contact plate, with bulb, without handle	D-001.78.106
XHL Xenon Halogen spare bulb	X-001.88.109



mini3000 Battery Handle with two dry batteries IEC LR6 size AA	D-001.79.021
---	---------------------

[HEINE mini 3000®]

	2.5 V
miniNT charger set for mini3000	
miniNT charger, 2 rechargeable batteries NiMH 2Z, 2 bottom inserts for mini3000 rechargeable handle	X-001.99.485
Rechargeable battery NiMH 2Z	X-001.99.487
miniNT adaptor set for mini2000	
2 bottom inserts for mini2000, 2 adaptors	X-001.99.491

Please note: To convert **mini 2000** handles to rechargeable, please order X-001.99.485 and X-001.99.491 together.

[HEINE mini 3000® SETS]

mini 3000 ENT Diagnostic Sets	
Set D-851 complete with mini3000 Fiber Optic (F.O.) Otoscope, 4 reusable tips, 10 AllSpec disposable tips in hard case with mini3000 battery handle with batteries	D-851.10.021
as above, but with mini3000 Otoscope	D-851.20.021
mini 3000 Ophthalmological Diagnostic Sets	
Set D-852 complete with mini3000 Ophthalmoscope in hard case, with mini3000 battery handle with batteries	D-852.10.021
mini 3000 Diagnostic Sets	
Set D-854 complete with mini3000 Ophthalmoscope, mini3000 Otoscope, mini3000 combi lamp, mini3000 straight Laryngeal Mirror, straight with plane and concave mirror, 1 set = 4 reusable tips, in hard case, with mini 3000 battery handle with batteries	D-854.10.021
as D-854, but with nose speculum, 10 mm Ø	D-853.10.021
Set D-873 complete with mini3000 Ophthalmoscope, mini3000 F.O. Otoscope, 1 set = 4 reusable tips, 5 each of 2.5 and 4mm Ø AllSpec disposable tips, in hard case, with 2 mini 3000 battery handles with batteries	D-873.11.021
as D-873, but with mini3000 Otoscope	D-873.21.021
Reusable tips (black) for mini 3000, mini 3000 F.O. Otoscopes	
1 set = 4 pcs. reusable tips SANALON S 2.4/3/4/5 mm dia.	B-000.11.111
Reusable tip 2.4 mm dia.	B-000.11.107
Reusable tip 3.0 mm dia.	B-000.11.108
Reusable tip 4.0 mm dia.	B-000.11.109
Reusable tip 5.0 mm dia.	B-000.11.110
Reusable nasal speculum 10 mm dia.	B-000.11.143
AllSpec Disposable tips (grey) for mini 3000 and mini 3000 F.O. Otoscopes	
Pack of 1000 2.5 mm dia. (infants) 1000 pcs./1 pack	B-000.11.128
Pack of 1000 4 mm dia. (adults) 1000 pcs./1 pack	B-000.11.127
Case of 10 dispenser packs of 250 tips, 2.5mm dia (infants)	B-000.11.151
Case of 10 dispenser packs of 250 tips, 4 mm dia (adults)	B-000.11.150
AllSpec Tip-dispenser	B-000.11.148



■ HEINE SUBSIDIARY COMPANIES
■ HEINE REPRESENTATIVES

A company makes history and gives a name to outstanding product quality: HEINE.

Founded in 1946, family owned, rich in tradition and today international market leader in a complete range of diagnostic instruments. HEINE products meet all international standards (ISO/CE) and represent the leading edge in precision and ergonomic design. Ongoing dedication to research and development has created a solid and diverse technology base for continued quality and product leadership. Over 500 employees worldwide contribute to this success.

www.heine.com

GERMANY

HEINE Optotechnik GmbH & Co. KG
82211 Herrsching
Kientalstr. 7
Tel. +49 (0) 81 52-38 0
Fax +49 (0) 81 52-3 82 02
E-Mail: info@heine.com

NORTH AMERICA

HEINE USA LTD.
10 Innovation Way
Dover, NH 03820
Tel. (603) 7 42-71 03
Fax (603) 7 42-72 17
Toll Free (800) 367- 4872
E-Mail: service@heine-na.com

AUSTRALIA

HEINE AUSTRALIA PTY. LTD.
Unit 9, 98 Old Pittwater Road
P.O. Box 717 Brookvale 2100
New South Wales
Tel. +61 (0) 2-99 38 95 00
Fax +61 (0) 2-99 39 23 05
E-Mail: info@heine.com.au

SWITZERLAND

HEINE (Schweiz) AG
Tobelackerstr. 9
CH-8212 Neuhausen am Rheinfall
Tel. +41 (0) 52-6 72 22 66
Fax +41 (0) 52-6 72 63 77
E-Mail: heine-schweiz@schaffhausen.ch

We reserve the right to change specification without notice.

Your distributor

10/09. A-000.00.235 e

