

RTA - 1800

Cleaning table with ultrasonic sink



► **Dimensions:** L 1800 x D 605 x H 900 (mm)

- Made entirely in stainless steel AISI 304, with double-wall doors containing fireproof sound-damping material and silicon edge seal that guarantees a perfect door seal and valid barrier to dust and dirt.
- Equipped of washing sink 500x400xh250mm for pre-cleaning of processed articles using a shower tap with cold and hot water between a demineralized water spray gun and a compressed air spray gun.
- Equipped of 20 Lt ultrasonic sink 510x300xh200mm configured at 35 KHz with the possibility to adjust the power and time, with manual drainage system.